

# Conometric retention and dynamic loading over time

Tebbel, Florian  
Halldin, Anders  
Frotscher, Marcel

## Abstract

Retention of the conometric connection over time has been investigated under dynamic loading in a test set-up based on ISO14801:2016.


The samples were subjected to dynamic load for  $5 \times 10^6$  cycles in total. Retention was measured after 1,  $0.35 \times 10^6$ ,  $1.25 \times 10^6$ ,  $2.5 \times 10^6$ ,  $3.75 \times 10^6$  and  $5 \times 10^6$  load cycles. The results showed that the ratio between axial load and retention increased over time but flattened out on a level between 0.6 – 0.85 after approx.  $1.25 \times 10^6$  load cycles.

It was concluded that the conometric connection of the tested samples maintains its stability over time under cyclic loading. There was no effect of wear detected that had a negative influence on the conometric connection.



Static & dynamic test setup

## Results



Retention during cyclic loading

The ratio increases over time but flattens out on a level between 0.6 – 0.85 after approx.  $1.25 \times 10^6$  cycles. No statistical difference of ratio between  $1.25 \times 10^6$  –  $5 \times 10^6$ , however, there is a statistically significant difference between cycle 1 and  $0.35 \times 10^6$  –  $5 \times 10^6$ . There is sufficient evidence to assume that the data is normally distributed (probability plot  $p \geq 0.05$ ).

## Background and Aim

An alternative to using cement- or screw-retained coupling for single tooth restorations on implants, is a conometric connection. This excludes the risk for excess cement or the disadvantage with a screw access hole. The retention of a single tooth friction-retained connection is dependent on the cone angle, coefficient of friction, initial push-in force and the external load situation. This study was set up to evaluate the conometric mode of retention over time under dynamic loading.

## Conclusion

It was concluded that the conometric connection of the tested samples maintains its stability over time under cyclic loading. There was no effect of wear detected that had a negative influence on the conometric connection.

## Methods and Materials

In modified  $30^\circ$  test setup, acc. to ISO 14801:2016, fully embedded implants, with conometric straight Acuris abutments and caps, were loaded with defined dynamic load at 15 Hz for  $5 \times 10^6$  cycles. Retention of the conical connection abutment to cap was tested after 1,  $0.35 \times 10^6$ ,  $1.25 \times 10^6$ ,  $2.5 \times 10^6$ ,  $3.75 \times 10^6$  and  $5 \times 10^6$  loading cycles in a pull-off test with a constant testing speed of 2 mm/min. The ratio retention/axial load was used to evaluate the retention over time under cyclic loading. The axial load was defined as follows: axial load = fatigue load  $\times \cos(30^\circ)$ .

## References