

Your Order Form for TiBase, Abutment Screw and ScanPost



Company name	
Surname, first name	
House number, street	
Zip code, City	
Phone	Fax
E-Mail	
Customer number	

TiBase and ScanPost are delivered without Scanbody. Please order this separately!



















The following components are compatible depending on the connection

Con- nection	Scanbodies for Omnica ⁴⁾		Scanbodies for Bluecam ⁴⁾		inCoris ZI meso F0.5		inCoris ZI meso F2	
	REF	Quantity in units	REF	Quantity in units	REF	Quantity in units	REF	Quantity in units
S	6431311		6431295		6231802		6231828	
L	6431329		6431303		6231810		6231836	


¹⁾ 1x Titanium base, 1x Abutment screw, ²⁾ 2x Abutment screw, ³⁾ 1x ScanPost, 1x Abutment screw, ⁴⁾ 36x Scanbody

Manufacturer / Implant	Implant diameter	Platform	Con- nection	TiBase ¹⁾			Abutment Screw ²⁾			ScanPost ³⁾		
				REF	Quantity in units	Tightening torque	REF	Quantity in units	REF	Quantity in units		
Dentsply Sirona Implants AstraTech Osseospeed EV	3		3.0	S	AT EV 3.0 GH1 S	6586304	25 Ncm	AT EV 3.0	6586262		AT EV 3.0 S	6586353
	3,6		3.6	S	AT EV 3.6 GH1 S	6586312		AT EV 3.6	6586270		AT EV 3.6 S	6586361
	4,2		4.2	L	AT EV 4.2 GH1 L	6586320		AT EV 4.2	6586288		AT EV 4.2 L	6586379
	4,8		4.8	L	AT EV 4.8 GH1 L	6586338		AT EV 4.8	6586296		AT EV 4.8 L	6586387
	5,4		5.4	L	AT EV 5.4 GH1 L	6586346		AT EV 5.4	6593714		AT EV 5.4 L	6586395
Dentsply Sirona Implants AstraTech Osseospeed TX	3.5 S / 4.0 S		3.5 / 4.0	L	AT OS 3.5/4.0 L	6282532	25 Ncm	AT OS 3.5/4.0	6460344		AT OS 3.5/4.0 L	6431055
	4.5 / 5.0 / 5.0 S		4.5 / 5.0	L	AT OS 4.5/5.0 L	6282540		AT OS 4.5/5.0	6460443		AT OS 4.5/5.0 L	6431063
Dentsply Sirona Implants Ankylos	A, B, C, D	C/X	S		ANK C/ GH1 S	6586528	15 Ncm	Not available		ANK S	6586569	
					ANK C/ GH2 S	6586536						
					ANK /X GH1 S	6586544						
					ANK /X GH2 S	6586551						
Dentsply Sirona Implants Frialit / Xive	3,4		3.4	S	FX 3.4 S	6282433	25 Ncm	FX 3.4, 3.8, 4.5, 5.5	6460476		FX 3.4 S	6430891
	3,8		3.8	S	FX 3.8 S	6282441					FX 3.8 S	6430909
	4,5		4.5	L	FX 4.5 L	6282458					FX 4.5 L	6430917
	5,5		5.5	L	FX 5.5 L	6282466					FX 5.5 L	6430925



Manufacturer / Implant	Implant diameter	Platform	Con-nection	TiBase ¹⁾				Abutment Screw ²⁾			ScanPost ³⁾		
				REF	Quantity in units	Tightening torque	REF	Quantity in units	REF	Quantity in units			
Biomet 3i Certain* (Internal connection)	3,4	 3.4	S	B C 3.4 S	6308048	20 Ncm	B C 3.4, 4.1, 5.0	6460450	B C 3.4 S	6431212			
	4,1	 4.1	L	B C 4.1 L	6308097				B C 4.1 L	6431220			
	5	 5.0	L	B C 5.0 L	6308121				B C 5.0 L	6431238			
Biomet 3i External hex	3,4	 3.4	L	B O 3.4 L	6282557	35 Ncm	B O 3.4, 4.1, 5.0	6460468	B O 3.4 L	6431089			
	4,1	 4.1		B O 4.1 L	6282565				B O 4.1 L	6431105			
	5	 5.0		B O 5.0 L	6282573				B O 5.0 L	6431113			
BioHorizons (Internal connection)	3,0 / 3,8	 3.0	S	BH 3.0 S	6532779	30 Ncm	BH 3.0	6561240	BH 3.0 S	6532761			
	3,0/3,5/3,8/4,0/4,6	 3.5	L	BH 3.5 L	6532894				BH 3.5 L	6532886			
	4,0/4,6/5,0/5,8	 4.5	L	BH 4.5 L	6532951				BH 4.5 L	6532944			
	5,0/5,8/6,0	 5.7	L	BH 5.7 L	6536242				BH 5.7 L	6536234			
Nobel Biocare Replace (3-channel internal connection)	3,5	 NP	L	NB RS 3.5 L	6282474	35 Ncm	NB RS 3.5	6460526	NB RS 3.5 L	6430933			
	4,3	 RP		NB RS 4.3 L	6282482				NB RS 4.3 L	6430941			
	5	 WP		NB RS 5.0 L	6282490				NB RS 4.3, 5.0, 6.0	6460534	NB RS 5.0 L	6430958	
	6	 6.0		NB RS 6.0 L	6282508				NB RS 6.0 L	6430982			
Nobel Biocare Nobel Active (Internal conical connection)	3,5	 NP	L	NB A 4.5 L	6308188	25 Ncm	NB A 4.5	6460484	NB A 4.5 L	6431279			
	4,3 / 5,0	 RP		NB A 5.0 L	6308253	35 Ncm	NB A 5.0	6460492	NB A 5.0 L	6431287			
Nobel Biocare Brånemark*	3,3	 NP	L	NB B 3.4 L	6282516	35 Ncm	NB B 3.4	6460500	NB B 3.4 L	6431006			
	3,75 / 4,0	 RP		NB B 4.1 L	6282524				NB B 4.1	6460518	NB B 4.1 L	6431022	



Manufacturer / Implant	Implant diameter	Platform	Con-nection	TiBase ¹⁾				Abutment Screw ²⁾				ScanPost ³⁾			
				REF	Quantity in units	Tightening torque	REF	Quantity in units	REF	Quantity in units					
Osstem (USA: Hiossen) Osstem TS (USA: Hiossen ET)	3,5	 Mini	L	O TS 3.5 L	6527035	20 Ncm	O TS 3.5	6561208		O TS 3.5 L	6534197				
	4,0/4,5/5,0/6,0/7,0	 Standard		O TS 4.0 L	6527043		30 Ncm	O TS 4.0		6561232	O TS 4.0 L		6536846		
Straumann Standard Tissue Level	3,3	NN (3,5 mm)	L	S SO 3.5 L	6284231	35 Ncm	S SO 3.5	6460559		S SO 3.5 L	6431162				
	3,3 / 4,1 / 4,8	RN (4,8 mm)		S SO 4.8 L	6284249		S SO 4.8, 6.5	6460567		S SO 4.8 L	6431170				
	4,8	WN (6,5 mm)		S SO 6.5 L	6284256		S SO 6.5 L	6431196							
Straumann Bone Level	3,3	NC (3,3 mm)	L	S BL 3.3 L	6308154	35 Ncm	S BL 3.3, 4.1	6460542		S BL 3.3 L	6431246				
	4,1 / 4,8	RC (4,1 mm / 4,8 mm)		S BL 4.1 L	6308337		S BL 4.1 L	6431253							
Thommen Medical SPI Element, SPI Contact, SPI Element Inicell, SPI Contact Inicell	3,5	 3,5	S	TM 3.5 S	6531854	25 Ncm	TM 3.5, 4.0, 4.5, 5.0, 6.0 6561273		TM 3.5 S	6544386					
	4	 4		TM 4 S	6532829				TM 4 S	6544394					
	4,5	 4,5		TM 4.5 S	6532837				TM 4.5 S	6544402					
	5	 5		TM 5 S	6544360				TM 5 S	6544410					
	6	 6		TM 6 S	6544378				TM 6 S	6544428					
Zimmer Tapered Screw-Vent	3,7 / 4,1	 3,5	L	Z TSV 3.5 L	6282581	30 Ncm	Z TSV 3.5, 4.5, 5.7 6460575		Z TSV 3.5 L	6431139					
	4,7	 4,5		Z TSV 4.5 L	6282599				Z TSV 4.5 L	6431147					
	6	 5,7		Z TSV 5.7 L	6282607				Z TSV 5.7 L	6431154					
Medentika M-Implant	3.5 / 4.0 / 4.5 / 5.0	3.5/5.0	L	MI 3.5 / 5.0 L	6308295	25 Ncm	Not available		Not available						