



Rinn®
Intraoral Imaging Accessories

Product Catalog 2020

THE DENTAL
SOLUTIONS
COMPANY™

 Dentsply
Sirona

Overview

Rinn® XCP-ORA® One Ring & Arm Positioning System.....	3	Film	
		Universal & Autoclavable Film Holders	
		Rinn® XCP®	11
		Rinn® Flip-Ray®	10
		Rinn® EndoRay® II	11
Digital Sensors		Radiation Protection	
Universal & Autoclavable Sensor Holders		Soothe-Guard Air®	12
Rinn® XCP-DS Fit®	4, 5	Lead-Free Dental X-Ray Aprons	12
Rinn® Snap-A-Ray® DS	6	Soothe-Guard®	13
Rinn® Snap-A-Ray® DS Endo	6	Lead-Lined Dental X-Ray Aprons	13
Rinn® Eeze-Grip®	6		
Universal & Disposable Sensor Holders		Rinn User Guide	15
Rinn® Uni-Grip® 360	7	Rinn Color Guide	15
Sensor Covers	5		
Phosphor Plates			
Universal & Autoclavable Phosphor Plate Holders			
Rinn® XCP-PSP Fit®	8		
Universal & Disposable Phosphor Plate Holders			
Rinn® Uni-Grip®	7		
Phosphor Plate Envelopes	9		

Decrease clutter with only **1 Arm and 1 Ring.**




Rinn® XCP-ORA® One Ring & Arm Positioning System

**One Ring. One Arm. One Goal:
Faster & Easier Radiographic Accuracy.**


- 3-in-1 positioning convenience for anterior, posterior, and bitewing
- Reduces components needed to accurately position radiographic media
- To be used with Rinn® biteblocks for
 - Digital Sensors (XCP-DS Fit®)
 - Phosphor Plates (XCP-PSP Fit®)
 - Film (XCP®)
- Works with round and rectangular collimators
- Autoclavable

ORA Ring & Arm	Qty
550771 XCP-ORA® Arm & Ring	1 pc.
550772 XCP-ORA® Arm & Ring	2 pcs.
550773 XCP-ORA® Ring	1 pc.
550774 XCP-ORA® Arm	1 pc.


For Digital Sensors:



For Phosphor Plates:



For Film:



Tips & Tricks on Sensor Positioning

Would you like to optimize your intraoral images while using digital sensors?

Check out our Youtube Channel and discover helpful tutorials around sensor positioning.


www.youtube.com

- 1
- 2
- 3



Please note: The green Endo-Biteblocks of the product lines Rinn® XCP®, XCP-PSP Fit™ and XCP-DS Fit® need the XCP® positioning arm and ring (see page 5, 9 and 11).

Universal & Fast Loading Design to fit all Digital Sensors.



Standard

Endo



Rinn® Snap-A-Ray® DS Universal Sensor Holder

Fast and easy to use.
For intraoral radiographs with bisecting technique.

- Provides accurate positioning
- Sliding lock holds 1, 1.5, and 2 securely
- For all anterior and posterior radiographs
- Autoclavable
- Snap-A-Ray® DS Endo features a cut out bite surface to allow room for your work




Use with the optional 2-in-1 Arm and Ring kit for alignment.

Rinn® Snap-A-Ray® DS Universal Sensor Holder	Qty
550321 Snap-A-Ray® DS	1 pc.
550323 Snap-A-Ray® DS	3 pcs.
550421 Snap-A-Ray® DS Endo	1 pc.
550423 Snap-A-Ray® DS Endo	3 pcs.

Rinn® Snap-A-Ray® DS Positioning Arm and Ring Kit	Qty
550322 Snap-A-Ray® DS Positioning Arm and Ring Kit	1 pc.
Contains each: 2-in-1 aiming ring (1), Straight arm (1), Angled arm (1)	

Rinn® Eezee-Grip® Digital Sensor Holder

- Fits most sensor brands and sizes
- Great for all anterior and posterior radiographs, including the third-molar region
- Patented flexible grips designed for patient comfort
- Autoclavable



Rinn® Eezee-Grip® Digital Sensor Holder	Qty
550293 Eezee-Grip® Digital Sensor Holder	1 pc.
550294 Eezee-Grip® Digital Sensor Holder	3 pcs.

Rinn® Uni-Grip® 360 Universal Sensor Holder

84% of sensor users found Uni-Grip® 360 to be faster than individual holders.¹

- Designed to capture an entire full mouth series using only one holder
- Single-use holder accommodates any #1, #1.5 or #2 sized digital sensor
- Easy rotating tab for posteriors, anteriors and bitewings
- Perforated adhesive tabs for fast, easy removal
- Enhanced patient comfort due to compact size and round corners



Rinn® Uni-Grip® 360	Qty
550052 Uni-Grip® 360 Kit	1 pc.
Contains: 50 Uni-Grip® 360 Universal Sensor Holders, 50 Sensor Covers	
550252 Uni-Grip® 360 Combo Pack	1 pc.
Contains: 100 Uni-Grip® 360 Universal Sensor Holders, 100 Sensor Covers, 1 Uni-GripAR® Positioning Arm and Ring	

Rinn® Uni-Grip® Universal Sensor Holder for Film and Phosphor Plates



Rinn® Uni-Grip®	Qty
550055 Uni-Grip® Kit for Film and Phosphor Plates	1 pc.
Contains: package of 50 Uni-Grip® Film and Phosphor Plate Holders	
550255 Uni-Grip® Combo Kit for Film and Phosphor Plates	1 pc.
Contains: two 550055 kits, one Uni-GripAR®	


Sliding the bitepiece along the channels enables you to take anterior, posterior and bitewing x-rays.



¹ Data on file.

Less Shifting. More Security.


Designed to keep phosphor plates in place.



Rinn® XCP-PSP Fit® Phosphor Plate Holder System


Less shifting and movement - minimizing the risk of errors and retakes

- Secure band envelopes work together with the biteblocks to keep the phosphor plate in place.
- Fits all major phosphor plate brands in sizes 0, 1, 2 and 3
- Protects phosphor plates
- Simplifies with color-coding assembly
- Enhances patient comfort with soft envelopes and rounded biteblocks
- Decreases clutter when used with Rinn® XCP-ORA® One Ring & Arm Positioning System



Rinn® XCP-PSP Fit® Phosphor Plate Envelopes with Secure Band

Envelopes with Secure Band		Qty
552714	Phosphor Plate Envelopes, Size #2	100 pcs.
552710	Phosphor Plate Envelopes, Size #0	100 pcs.
552711	Phosphor Plate Envelopes, Size #1	100 pcs.
552712	Phosphor Plate Envelopes, Size #2	200 pcs.
552713	Phosphor Plate Envelopes, Size #3	100 pcs.



Rinn® XCP-PSP Fit® Biteblocks



Anterior Thick		Fits Phosphor Plate Size		Qty	Bitewing Horizontal		Fits Phosphor Plate Size		Qty
552730	Anterior Thick	1		2 pcs.	552760	Bitewing Horizontal	0		4 pcs.
552731	Anterior Thick	2		2 pcs.	552761	Bitewing Horizontal	1		4 pcs.
Anterior Thin					552762	Bitewing Horizontal	2		4 pcs.
552740	Anterior Thin	0		4 pcs.	552763	Bitewing Horizontal	2		10 pcs.
552741	Anterior Thin	1		4 pcs.	552764	Bitewing Horizontal	3		4 pcs.
552742	Anterior Thin	2		4 pcs.	Bitewing Vertical				
552743	Anterior Thin	2		10 pcs.	552770	Bitewing Vertical	2		4 pcs.
Posterior					Endodontic				
552750	Posterior	0		4 pcs.	552781	Endo Horizontal	2		4 pcs.
552751	Posterior	1		4 pcs.	552790	Endo Vertical	1		4 pcs.
552752	Posterior	2		4 pcs.	552791	Endo Vertical	2		4 pcs.
552753	Posterior	2		10 pcs.					

Rinn® XCP-ORA® One Ring & Arm Positioning System



Rinn® XCP® Aiming Ring & Arm



Rinn® XCP-ORA® Ring & Arm		Qty	Rinn® XCP® Endo Arm & Ring		Qty
550771	XCP-ORA® Arm & Ring	1 pc.	550597	Endo Arm (No pins)	1 pc.
550772	XCP-ORA® Arm & Ring	2 pcs.	550598	Endo Ring (Green)	1 pc.
550773	XCP-ORA® Ring	1 pc.			
550774	XCP-ORA® Arm	1 pc.			

Rinn® XCP-PSP Fit® Kits Phosphor Plate Holder System




Rinn® XCP-PSP Fit® Phosphor Plate Holder System for Endo



Rinn® XCP-PSP Fit® Kits		Rinn® XCP-PSP Fit™	
552700	Complete Size 0 Solution Kit Contains: Envelopes #0 (100), Anterior Thin #0 (4), Posterior #0 (4), Bitewing Hor. #0 (4), ORA Ring & Arm (1)	552703	Complete Endo Kit Contains: Envelopes #1 (100), Envelopes #2 (100), Endo Vert. #2 (4), Endo Hor. #2 (4), Endo Vert. #1 (4), XCP Endo Arm (1), XCP Endo Ring (1)
552701	Complete Size 1 Solution Kit Contains: Envelopes #1 (100), Anterior Thick #1 (2), Anterior Thin #1 (4), Posterior #1 (4), Bitewing Hor. #1 (4), ORA Ring & Arm (1)		
552702	Complete Size 2 Solution Kit Contains: Envelopes #2 (100), Anterior Thick #2 (2), Anterior Thin #2 (8), Posterior #2 (8), Bitewing Hor. #2 (8), ORA Ring & Arm (1)		

Angled for
hard-to-reach areas



Rinn® Snap-A-Ray® Xtra Film & Phosphor Plate Holder

Film and Phosphor Plate Holder for Bisecting Technique

- Cushioned media grips hold film and phosphor plates securel
- Cushioned bite areas designed to enhance patient comfort
- Media can be angled for hard-to-reach areas
- Autoclavable

Rinn® Snap-A-Ray® Xtra Film & Phosphor Plate Holder	Qty
540501 Snap-A-Ray® Xtra Film & Phosphor Plate Holder	1 pc.
540503 Snap-A-Ray® Xtra Film & Phosphor Plate Holder	3 pcs.

Rinn® Flip-Ray® Universal Film & Phosphor Plate Holder

Flip media holder from horizontal to vertical!

- Adjustable ring and media holder flip orientations to change quadrants
- Accommodates size 0, 1, and 2 film and phosphor plates
- Available with adult size or smaller, pedo size media holder
- Autoclavable

Rinn® Flip-Ray®	Qty
552600 Flip-Ray® Complete Kit: 2 Flip-Ray PA (1 Adult, 1 Pedo) and 2 Flip-Ray BW	4 pcs.
552601 Flip-Ray® PA, Pedo	4 pcs.
552602 Flip-Ray® PA, Adult	4 pcs.
552603 Flip-Ray® BW Universal Bitewing Holder	5 pcs.




Rinn® Flip-Ray® PA
for Anterior and Posterior



Rinn® Flip-Ray® BW
for Bitewings

Rinn® XCP® Holder System

Recommended for Film



Rinn® XCP® Kits	
542001 XCP® Evolution 2000 Kit Anterior, Posterior and Bitewings	Contains: Anterior (4), Posterior (8), Bitewing Size #0 & #1, hor. (1), Bitewing Size #1, vert. (1), Bitewings Size #2, hor. (2), Bitewings Size #2, vert. (1), Positioning Arms & Rings: Anterior (1), Posterior (2), Bitewing (1), Positioning Poster (1)
542003 XCP® Evolution 2000 Kit Complete Kit (incl. Endo)	Contains: all components of 542001 as well as Positioning Arm Endo (1), Positioning Ring Endo (1), Biteblocks Endo (4)
542045 XCP® Endo Kit	Contains: XCP® Positioning Arm Endo (1), XCP® Positioning Ring Endo (1), XCP® Biteblocks Endo (3)

XCP® Biteblocks	Qty
540861 Anterior	25 pcs.
540862 Posterior	25 pcs.
540867 Anterior	6 pcs.
540868 Posterior	6 pcs.
540863 Anterior and Posterior	10 pcs.
540929 Bitewings Size #2, Horizontal	6 pcs.
540931 Bitewings Size #3, Horizontal	6 pcs.
540935 Endo	3 pcs.

XCP® Positioning Arms	Qty
540857 Anterior	1 pc.
540858 Posterior	1 pc.
540927 Bitewings	1 pc.
550597 Endo (No pins)	1 pc.

XCP® Positioning Rings	Qty
540865 Anterior	1 pc.
540860 Posterior	1 pc.
540934 Bitewings	1 pc.
550598 Endo	1 pc.

EndoRay® II Film Holder

Fits over endodontic files, clamps and rubber dams.


- Use for endodontic radiography without the cheek interfering with the aiming arm
- Use with #1 or #2 film and phosphor plates
- Autoclavable

EndoRay® II	Qty
540303 EndoRay® II	2 pcs.
540307 EndoRay® II	4 pcs.



Rinn® XCP-ORA® Complete Kit for Film

Rinn® XCP-ORA®	
542071 XCP-ORA® Complete Kit for Film (no endo included)	Contains: 72 Biteblocks - Anterior (20) Posterior (20) Bitewings Size #0 & #1, Horizontal (10) Bitewings Size #1, Vertical (6) Bitewings Size #2, Horizontal (10) Bitewings Size #2, Vertical (6) ORA® Ring & Arm (1) Positioning Poster (1)





*Please note that due to a variety of monitor settings, final print and or other marketing collateral may slightly vary from the actual product color.

Introducing the improved

Soothe-Guard Air® Lead-Free Dental X-Ray Aprons

Soothe-Guard® Lead-Lined Dental X-Ray Aprons

Protecting your patients from unnecessary radiation exposure is a key clinical responsibility. Keeping them comfortable is key to ensuring their satisfaction and bringing them back to your office. The improved Rinn® Dental Aprons deliver on those expectations in a premium innovative design.

Innovative Design and Magnetic Closure Option

- Up to 10 styles for intraoral & extraoral imaging
- 2 lead equivalencies to protect your patients: Lead-Free or Lead-Lined - 0.35 mm and 0.5 mm
- 6 premium colors to complement your practice decor
- Innovative magnetic closure provides a secure fit and prevents hair and clothing from being entangled in the closure.
- Complies with global standards for radiation protection (IEC 61331-1:2014, PPE 2016/425).



Soothe-Guard Air® Lead-Free Aprons

Soothe-Guard Air Lead-Free aprons are up to 15% lighter than Soothe-Guard Lead-Lined aprons.

.35MM Lead Equivalency	Closure Type	Weight (lbs.)	Navy Blue	Light Blue	Cool Gray	Sand	Royal Blue	Charcoal
A Adult Apron with Collar	Magnetic	5.95	861805001	861905001	862005001	862105001	862205001	862305001
B Child Apron with Collar	Magnetic	3.31	861804701	861904701		862104701	862204701	
C Adult Apron Ext. Collar	Magnetic	6.17				862106001		862306001
D Child Apron Ext. Collar	Magnetic	3.53				862105701		862305701
E Adult Apron without Collar	No closure	5.18	861804902	861904902	862004902	862104902	862204902	862304902
F Thyroid Collar	Magnetic	0.79						862305101
G Adult Dual-Pano	Standard	9.15	861805400	861905400	862005400	862105400	862205400	862305400
H Child Dual-Pano	Standard	5.18	861805500	861905500		862105500	862205500	
.5MM Lead Equivalency								
A Adult Apron with Collar	Magnetic	7.61	871805001					

Soothe-Guard® Lead-Lined Aprons


















.35MM Lead Equivalency	Closure Type	Weight (lbs.)	Navy Blue	Light Blue	Cool Gray	Sand	Royal Blue	Charcoal
A Adult Apron with Collar	Standard	6.50	661805000	661905000	662005000	662105000	662205000	662305000
B Child Apron with Collar	Standard	3.97	661804700	661904700		662104700	662204700	
C Adult Apron Ext. Collar	Standard	6.73				662106000		662306000
E Adult Apron without Collar	No closure	5.84	661804902	661904902	662004902	662104902	662204902	662304902
.5MM Lead Equivalency								
E Adult Apron without Collar	No closure	7.59	671804902					





According to the ADA:

“ All protective shields should be evaluated for damage (ie. tears, folds, and cracks) monthly using visual and manual inspection.”²

Note: A crack in the outer layer could expose the lead inside, which could transfer to the patient. We recommend **replacing your aprons every 3 years.**

² Data on file.

<p>Media in Use</p>	<p>Digital Sensors</p>  <p>e.g. Dentsply Sirona Xios</p>	<p>Phosphor Plates</p>  <p>e.g. Dentsply Sirona Xios</p>	<p>Film</p>  <p>e.g. Carestream E-Speed</p>
<p>General Use</p>	 <p>Rinn® XCP-DS Fit™ Page 4/5</p>	 <p>Rinn® XCP-PSP Fit™ Page 8/9</p>	 <p>Rinn® XCP® Page 11</p>
<p>Endodontic Procedures</p>	 <p>Rinn® XCP-DS Fit® Endo Kit Page 4/5</p>  <p>Rinn® Snap-A-Ray® DS Endo Page 6</p>	 <p>Rinn® XCP-PSP Fit™ Endo Kit Page 8/9</p>	 <p>Rinn® XCP® Endo Kit Page 11</p>
<p>Pediatric Procedures</p>	 <p>Rinn® XCP-DS Fit® Size 1 Vertical Bitewing to be used with size 0 sensors for horizontal Bitewings Page 4/5</p>	 <p>Rinn® XCP-PSP Fit™ Size 0 & Size 1 Kits Page 8/9</p>	 <p>Rinn® Flip-Ray® Pedo Page 11</p>
<p>Single Use Option</p>	 <p>Rinn® Uni-Grip® 360 Page 7</p>	 <p>Rinn® Uni-Grip® Page 7</p>	
<p>Bisecting Options (for challenging situations such as Narrow Jawbone and Tori)</p>	 <p>Rinn® Snap-A-Ray® DS Page 6</p>	 <p>Rinn® Snap-A-Ray® Xtra Page 10</p>	

<p>Anterior</p>		
<p>Posterior</p>		
<p>Bitewing</p>		
<p>Endo</p>		

For more information visit www.dentsplysirona.com.

Find us on  www.facebook.com/dentsplysirona.preventive

© 2020 Dentsply Sirona
RIN08-0619-1-EU Rev. 1 (0220)

THE DENTAL
SOLUTIONS
COMPANY™

