

Learn more about our resto workflow and other workflows at Dentsply Sirona: dentsplysirona.com/clinictolab

Digital restorative workflow with Primescan



Better everyday restorative dentistry

Experience faster, easier, more consistent high-quality outcome in all restorative procedures thanks to Primescan powered by DS Core™. Digitize your existing restorative workflows or expand to new restorative treatment offers from single to multiple units, including large aesthetics and more complex cases.

Digital Impression with Primescan

Outstanding quality

- Highly accurate full arch scan in less than a minute
- Easy scanning of shiny materials and edentulous areas
- Scan as the foundation for better function, and esthetic

Comfort

- Enhance treatment experience and discuss cases in real time while the patient is still in the dental chair
- Intuitive use allows for an easy delegation to your team

Open Scanner

- Validated connection to all major lab CAD softwares for easy collaboration




Freedom of collaboration with Connect Case Center

- With one click, your lab of choice can access to the scan data via Connect Case Center Inbox.
- Your partner can easily proceed with the design regardless of the CAD software used for a shorter treatment time
- Anonymized and immediate data transfer within seconds (HIPAA and GDPR compliance)
- No license fees



Benefits of DS Core™ for Digital Impression

With DS Core™, Dentsply Sirona offers a cloud-based platform that combines products, services and technology that empowers your practice growth through efficiency and connectivity.



Share & collaborate

Efficient collaboration with specialists, partners and labs. Send patient records to your colleagues with the click of a button.

One central library

Easily combine digital impressions and x-rays for any workflow in one library. Scans and images will be automatically uploaded from your Dentsply Sirona equipment.

Easy access

Access the platform from multiple devices easily through a web browser. View, download and manage case files from different equipments, no matter where they are located.

The restorative workflow step-by-step



Scanning

Primescan can quickly and accurately take a full-arch scan. Digital impression has the advantage that you can interrupt and resume the scan procedure at any time. Additionally, you can assess the quality of the model right there on the screen. This procedure can be delegated, saving time and bringing more joy to your practice team.

Using
DS CORE™

A copy of the model data is automatically uploaded and added to your patient records as a .dxd file. You can access these files remotely using a web browser on any scanner or workstation at any time.*



Data transmission

Send the scan data via the Connect Case Center to any lab for easy and safe ordering with one click. Additionally, you can make use of DS Core™ to share case files with colleagues and to collaborate on treatment plans. You can also send all patient records from one platform in case of referring your patient to a specialist.


Using
DS CORE™

Experience a new way of collaboration with peers or specialists within digital workflows and complex restorative cases: Easily combine digital impressions, x-rays, patient videos or photos all from one library.*



Data reception by the lab of your choice

All third parties are able to access scan data transmitted via Connect Case Center by installing a Connect Case Center Inbox. You and your lab benefit from validated connections to all major lab CAD softwares. Your collaboration will become closer than ever before.



Treating the patient

With Primescan you can serve a broad variety of restorative indications giving your patients the treatment they deserve. Provide them with indirect restorations such as crowns, inlays, onlays and bridges, as well as implant restorations such as abutments and screw-retained crowns. Additionally, you can include esthetic treatment options like veneers and smile makeovers to your practice offer.

“This evolution of partnership is changing. It’s stronger and like so many things that are online, it’s upon the practitioner’s convenience and their preferences and we’re here to assist them with whatever they desire.”

Devaughn Fraser | Highland Dental Laboratory | Calgary, Canada



* This is not a medical device software. It is not intended to diagnose, treat, cure, or prevent any disease or health condition. Use only for information, education or sharing purposes. Service is subject to country availability. Please contact with your local DS representatives.